BioPro Horizon Subtalar Implant

Indications for use:

The BioPro Horizon Subtalar Implant is indicated for use in the treatment of the hyperpronated foot and stabilization of the subtalar joint. It is designed to block forward, downward and medial displacement of the talus, thus allowing normal subtalar joint motion but blocking excessive pronation and the resulting sequela.

Conditions include:

- Flat foot treatment in children and adolescents
- Congenital flat foot
- Unsuccessful long term orthopedic treatment (orthotics)
- Tarsal coalitions
- Painfully flat foot
- Supple deformity in posterior tibial tendon dysfunction
- Paralytic flat foot
- Subtalar instability

Contra-indications:

1. A debilitating general health problem that might pose a significant threat to the life of the patient if subjected to a major surgical procedure
2. A previously infected joint that has not been quiescent for at least six months
3. A local or systemic infection (i.e. osteomyelitis)
4. Insufficient bone stock to support the prosthesis

The age of the patient must be balanced against the severity of the disability and the need for surgery

Sterile:

This is a single use device. Re-use of this device can result in the transfer of materials not limited to bone, tissue, blood or infectious disease. This device is provided sterile and re-sterilization of the device has not been validated.

Sterilized with ethylene oxide gas. *Caution: For one procedure only. Do not re-sterilize. Do not use if package is open or damaged.

The devices have been evaluated for safety and compatibility in the MR environment and are MR conditional. Contact BioPro for MR parameters. Contact surgeon if a change in performance or pain level is noticed

Instructions for use:

Step One:

Attention is directed to the sinus tarsi. A linear incision centering over the sinus tarsi is made into Relaxed Skin Tension Line (RSTL), approximately 2cm in length. The intermediate dorsal cutaneous nerve is identified, carefully retracted and preserved. Next, a linear incision is made into the retinaculum to expose the sinus tarsi.

Step Two:

A probe is used to slightly enlarge the sinus as well as the canalis tarsi and also to establish subtalar joint axis. The leading edge of the probe should be palpated at the medial aspect of the subtalar joint with slight tenting of the skin. This should be appreciated just inferior to the posterior tibial tendon and slightly inferior and anterior to the medial malleolus. Care is taken to preserve the interosseous talocalcaneal ligament.

Step Three:

Pass one of the provided guide wires through the cannulation on the Starter Probe until palpable under the skin on the medial aspect of the foot. Create a small incision over the point of the guide wire and pass through approximately 5mm of guide wire. Clamp a hemostat over the exposed guide wire to stabilize the wire. The Starter Probe can now be removed laterally.

Step Four:

Place the Sizer probes over the guide wire and insert into the Subtalar joint until the end of the probe is just past the longitudinal bisection of the talus, continue inserting the sequentially sized Sizer Probes until the desired restricted Subtalar joint motion and clinical correction is achieved. This can be assessed intraoperatively by evert and inverting the calcaneus and at the same time loading the lateral column. Approximately 4° to 6° of eversion of the calcaneus should be noted. Note the position of the skin line along the graduated markings on the sizing probe. Now, the Sizer Probe is removed and the guide wire is maintained.

Step Five:

Note: color-coded trial implants are available in the instrument kit, but it is recommended to progress direction into the final implant once proper correction is achieved with the Sizer Probes.

The appropriately sized final implant is chosen, inserted over the guide wire, and then threaded into place with the cannulated screwdriver (fig. 5). At this time the graduated markings of the screwdriver should match and correlate to the Sizer Probe placement markings.

Important: the sizer markings should not supersede the interoperative Subtalar joint motion evaluation of the surgeon.

Subtalar joint motion is evaluated and clinical correction is appreciated. At this time the placement of the implant can be appreciated (surgeon's discretion) with the C-arm on anterior to posterior ankle review.

Important: the positioning should be noted where the medial or leading edge of the implant should go slightly past the longitudinal bisection of the talus (1-2mm).

For a complete glossary of symbols, please visit www.bioproimplants.com/IFU. To receive a printed version within 5 business days, please contact orders@bioproimplants.com

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BioPro, Inc.