BioPro Toe M-P Joint – CE Toes

Indications for use:

- 1) Arthritic degradation of the metatarso-phalangeal joint of the great toe that has resulted in disabling pain, limited motion and loss of the normal ambulatory function of the forefoot.
- 2) Degenerative arthritis (hallux rigidus)
- 3) Rheumatoid Arthritis
- 4) Bunion deformity associated with arthritis of the great metatarso-phalangeal joint

Contra-indications:

- 1) A general health problem that might pose a significant threat to the life of the patient if subjected to a major surgical procedure.
- 2) An active infection or a previous infection of the lower extremity that has not been quiescent for at least six months.
- 3) A local or systemic infection.
- 4) Significant deficiency in the vascular supply to the extremity.
- 5) Severe structural deficiency of the sub-chondral bone that may result in insufficient support for the prosthesis.
- 6) A condition of the toe which may lend itself to a more conservative procedure.
- 7) Severe compromise of the supporting muscles or ligaments about the toe.

Sterile:

For a complete glossary of symbols, please visit www.bioproimplants.com/IFU. To receive a printed version within 5 business days, please contact orders@bioproimplants.com

Sterilized with ethylene oxide gas. Caution: For one procedure only. Do not re-sterilize. Do not use if package is open or damaged. This is a single use device. Re-use of this device can result in the transfer of materials not limited to bone, tissue, blood, or infectious disease. The device is provided sterile and re-sterilization of the device has not been validated. The devices have been evaluated for safety and compatibility in the MR environment and are MR conditional. Contact BioPro for MR parameters. Contact surgeon if a change in performance or pain level is noticed

Instructions for use:

Extreme care is to be taken to preserve all the capsular structures and to maintain attachment of the intrinsic musculature. The joint capsule is released subperiosteally from the base of the proximal phalanx at the medial, plantar and lateral aspects of the bone.

Using a narrow oscillating blade remove a section of bone sufficient to avoid prosthetic over spacing and **excessive joint tension** and to accommodate the thickness of the articulating plate of the implant.

Next, the intramedullary canal is prepared with a punch, the implant is sized with the most appropriately approximated dimension of the osteotomized phalanx, which does not extend beyond the margins, and is inserted into the remaining base of the proximal phalanx making sure that the implant is completely seated. Extra care is taken not to impinge any capsular structures.

At this time the metatarsal head is remodeled and all marginal osteophytes are completely resected. Range of motion of the joint should be smooth and without evidence of crepitation or binding.

Postoperative management:

The patient is allowed to ambulate with weightbearing to tolerance on the operated foot within limits imposed by postoperative discomfort, utilizing modified footwear. The progression to normal ambulation and the use of standard footwear is limited only by the persistence of postoperative swelling and discomfort.

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BioPro, Inc.