



Tray Layout Guide

The BioPro Foot Plating System tray utilizes a modular design allowing the user to customize the tray, thereby reducing the footprint in the operating room.

- Base Instrument Tray [ref #22925]
- 3.0/3.5 Plate Module [ref #22931]
- 2.8 Plate Module [ref #22930]
- Medial Column Plate Module [ref #22932]
- Calcaneal Plate Module [ref #22933]
- MTP Plate Module [ref #22929/23179]
- 2.0 Screw + Plate Module [ref #22926]
- 2.8 Screw Module [ref #22927]
- 3.0/3.5 Screw Module [ref #22928]
- Osteotomy Plate Module [ref #22934]




The Foot Plating System Main Tray [ref# 22925] is constructed with 3 levels. Each module is to be ordered separately.




Top Level

- 1 Plate Modules (ordered separately)
 - Medial Column Plate Module [ref #22932]
 - Calcaneal Plate Module [ref #22933]
- 2 Plate Module (ordered separately)
 - 3.0/3.5 Plate Module [ref #22931]
 - 2.8 Plate Module [ref #22930]
- 3 Drivers
- 4 Depth Gauges



Middle Level


- 5 3.0/3.5 Screw Module [ref #22928]
(ordered separately)
- 6 2.8 Screw Module [ref #22927]
(ordered separately)
- 7 2.0 Screw & Plate Module [ref #22926]
(ordered separately)
- 8 MTP Module [ref #22929]
(ordered separately)



Bottom Level

- 9 Osteotomy Plate Module [ref #22934]
(ordered separately)
- 10 Bending Irons
- 11 Double Drill Guides
- 12 Hohman Retractors
- 13 Miscellaneous Instruments (ordered separately)


2.8 Plate Module [ref# 22930]



2.8 Plate Module - 22930

Location	Screw Color	Item #	Description	Qty
1	•	20614	Straight Fracture Plate, 2.5-2.8mm, 4 Hole	1
2	•	20615	Straight Fracture Plate, 2.5-2.8mm, 5 Hole	1
3	•	20616	Straight Fracture Plate, 2.5-2.8mm, 6 Hole	1
4	•	20617	Straight Fracture Plate, 2.5-2.8mm, 7 Hole	1
5	•	20624	T-Shape Fracture Plate, 2.5-2.8mm, 3 Hole	1
6	•	20625	T-Shape Fracture Plate, 2.5-2.8mm, 4 Hole	1
7	•	20622	L-Shape Fracture Plate, 2.5-2.8mm, 3 Hole, Right	1
8	•	20623	L-Shape Fracture Plate, 2.5-2.8mm, 4 Hole, Right	1
9	•	20620	L-Shape Fracture Plate, 2.5-2.8mm, 3 Hole, Left	1
10	•	20621	L-Shape Fracture Plate, 2.5-2.8mm, 4 Hole, Left	1
11	•	20756	Medial TMT-1 Fusion Plate, 2.8mm, 4 Hole	1
12	•	20618	Y-Shape Fracture Plate, 2.5-2.8mm, 3 Hole	1
13	•	20619	Y-Shape Fracture Plate, 2.5-2.8mm, 4 Hole	1
14	•	20626	Cloverleaf Plate, 2.5-2.8mm, 3 Hole	1
15	•	20627	Cloverleaf Plate, 2.5-2.8mm, 4 Hole	1
16	•	20628	Open Wedge Plate, 2.5-2.8mm, 0mm Wedge	1
17	•	20629	Open Wedge Plate, 2.5-2.8mm, 2mm Wedge	1
18	•	20630	Open Wedge Plate, 2.5-2.8mm, 3mm Wedge	1
19	•	20631	Open Wedge Plate, 2.5-2.8mm, 4mm Wedge	1
20	•	20632	Open Wedge Plate, 2.5-2.8mm, 5mm Wedge	1
21	•	20633	Open Wedge Plate, 2.5-2.8mm, 6mm Wedge	1


3.0-3.5 Plate Module [ref# 22931]



3.0-3.5 Plate Module - 22931

Location	Screw Color	Item #	Description	Qty
1	● ●	20868	Straight Fusion Plate, 3.0-3.5mm, 2 Hole	1
2	● ●	20869	Straight Fusion Plate, 3.0-3.5mm, 4 Hole	1
3	● ●	20870	Straight Fusion Plate, 3.0-3.5mm, 6 Hole	1
4	● ●	21209	T-Shape Fracture Plate, 3.0-3.5mm, 3 Hole	1
5	● ●	20877	T-Shape Fracture Plate, 3.0-3.5mm, 4 Hole	1
6	● ●	20878	T-Shape Fracture Plate, 3.0-3.5mm, 6 Hole	1
7	● ●	21197	L-Shape Fracture Plate, 3.0-3.5mm, 3 Hole, Left	1
8	● ●	20871	L-Shape Fracture Plate, 3.0-3.5mm, 4 Hole, Left	1
9	● ●	20872	L-Shape Fracture Plate, 3.0-3.5mm, 6 Hole, Left	1
10	● ●	21199	L-Shape Fracture Plate, 3.0-3.5mm, 3 Hole, Right	1
11	● ●	20873	L-Shape Fracture Plate, 3.0-3.5mm, 4 Hole, Right	1
12	● ●	20874	L-Shape Fracture Plate, 3.0-3.5mm, 6 Hole, Right	1
13	● ●	21201	Cloverleaf Plate, 3.0-3.5mm, 3 Hole	1
14	● ●	20875	Cloverleaf Plate, 3.0-3.5mm, 4 Hole	1
15	● ●	20876	Cloverleaf Plate, 3.0-3.5mm, 6 Hole	1
16	● ●	20879	X-Plate, 3.0-4.0mm, Small	1
17	● ●	20880	X-Plate, 3.0-4.0mm, Medium	1
18	● ●	20881	X-Plate, 3.0-4.0mm, Large	1
19	● ●	20885	Rectangular Plate, 3.0-3.5mm, Small	1
20	● ●	20886	Rectangular Plate, 3.0-3.5mm, Medium	1
21	● ●	20887	Rectangular Plate, 3.0-3.5mm, Large	1

MTP Plate Module [ref# 22929/23179]




MTP Plate Module - 22929

Location	Screw Color	Item #	Description	Qty
1	• • •	20611	MTP Fusion Plate, 2.5-3.5mm, Right	1
2	• • •	20610	MTP Fusion Plate, 2.5-3.5mm, Left	1
3	• • •	20613	MTP Fusion Rev Plate, 2.8-3.5mm, Right	1
4	• • •	20612	MTP Fusion Rev Plate, 2.8-3.5mm, Left	1
5		20671	MTP Fusion Reamer, Convex, 16mm	1
6		20672	MTP Fusion Reamer, Convex, 18mm	1
7		20673	MTP Fusion Reamer, Convex, 21mm	1
8		20674	MTP Fusion Reamer, Convex, 24mm	1
9		20677	MTP Fusion Reamer, Concave, 16mm	1
10		20678	MTP Fusion Reamer, Concave, 18mm	1
11		20679	MTP Fusion Reamer, Concave, 21mm	1
12		20680	MTP Fusion Reamer, Concave, 24mm	1
13		17232	K-Wire .062" x 4" Non-Sterile	4

MTP Plate Module Gen 2 - 23179

Location	Screw Color	Item #	Description	Qty
1	• • •	23086	MTP Fusion Plate Narrow 5, 2.5-3.5mm, Right	1
2	• • •	23085	MTP Fusion Plate Narrow 5, 2.5-3.5mm, Left	1
1A	• • •	23084	MTP Fusion Plate Narrow 10, 2.5-3.5mm, Right	1
2A	• • •	20383	MTP Fusion Plate Narrow 10, 2.5-3.5mm, Left	1
3	• • •	20613	MTP Fusion Rev Plate, 2.8-3.5mm, Right	1
4	• • •	20612	MTP Fusion Rev Plate, 2.8-3.5mm, Left	1
5		20671	MTP Fusion Reamer, Convex, 16mm	1
6		20672	MTP Fusion Reamer, Convex, 18mm	1
7		20673	MTP Fusion Reamer, Convex, 21mm	1
8		20674	MTP Fusion Reamer, Convex, 24mm	1
9		20677	MTP Fusion Reamer, Concave, 16mm	1
10		20678	MTP Fusion Reamer, Concave, 18mm	1
11		20679	MTP Fusion Reamer, Concave, 21mm	1
12		20680	MTP Fusion Reamer, Concave, 24mm	1
13		17232	K-Wire .062" x 4" Non-Sterile	4


Osteotomy Plate Module [ref# 22934]



Osteotomy Module - 22934


Location	Screw Color	Item #	Description	Qty
1	● ●	20888	Evans Plate, 3.0-3.5mm, 4mm Step, Left	1
2	● ●	20889	Evans Plate, 3.0-3.5mm, 6mm Step, Left	1
3	● ●	20890	Evans Plate, 3.0-3.5mm, 8mm Step, Left	1
4	● ●	20891	Evans Plate, 3.0-3.5mm, 10mm Step, Left	1
5	● ●	20892	Evans Plate, 3.0-3.5mm, 4mm Step, Right	1
6	● ●	20893	Evans Plate, 3.0-3.5mm, 6mm Step, Right	1
7	● ●	20894	Evans Plate, 3.0-3.5mm, 8mm Step, Right	1
8	● ●	20895	Evans Plate, 3.0-3.5mm, 10mm Step, Right	1
9	● ●	20896	Dwyer Plate, 3.0-3.5mm, 4mm Step	1
10	● ●	20897	Dwyer Plate, 3.0-3.5mm, 6mm Step	1
11	● ●	20898	Dwyer Plate, 3.0-3.5mm, 8mm Step	1
12	● ●	20899	Dwyer Plate, 3.0-3.5mm, 10mm Step	1

Medial Column Plate Module [ref# 22932]




Medial Column Module - 22932				
Location	Screw Color	Item #	Description	Qty
1	● ●	20909	Medial Distal Column Plate, 3.0-3.5mm, Short	1
2	● ●	20908	Medial Distal Column Plate, 3.0-3.5mm, Long	1
3	● ●	20905	Medial 3 Column Plate, 3.0-3.5mm, Short, Left	1
4	● ●	20904	Medial 3 Column Plate, 3.0-3.5mm, Short, Right	1
5	● ●	20907	Medial 3 Column Plate, 3.0-3.5mm, Long, Left	1
6	● ●	20906	Medial 3 Column Plate, 3.0-3.5mm, Long, Right	1
7	● ●	20901	Medial 4 Column Plate, 3.0-3.5mm, Short, Left	1
8	● ●	20900	Medial 4 Column Plate, 3.0-3.5mm, Short, Right	1
9	● ●	20903	Medial 4 Column Plate, 3.0-3.5mm, Long, Left	1
10	● ●	20902	Medial 4 Column Plate, 3.0-3.5mm, Long, Right	1

Calcaneal Plate Module [ref# 22933]



Calcaneal Plate Module - 22933				
Location	Screw Color	Item #	Description	Qty
1	• •	20910	ORIF Calcaneal Plate, 3.0-3.5mm, Short	1
2	• •	20911	ORIF Calcaneal Plate, 3.0-3.5mm, Medium	1
3	• •	20912	ORIF Calcaneal Plate, 3.0-3.5mm, Long	1

2.0 Screw & Plate Module [ref# 22926]



2.0 Screw-Plate Module - 22926

Location	Screw Color	Item #	Description	Qty
1	●	21207	Straight Fracture Plate, 2.0mm, 3 Hole	1
2	●	20602	Straight Fracture Plate, 2.0mm, 4 Hole	1
3	●	20604	Straight Fracture Plate, 2.0mm, 6 Hole	1
4	●	20605	Straight Fracture Plate, 2.0mm, 8 Hole	1
5	●	21205	T-Shape Fracture Plate, 2.0mm, 3 Hole	1
6	●	20608	T-Shape Fracture Plate, 2.0mm, 4 Hole	1
7	●	20609	T-Shape Fracture Plate, 2.0mm, 6 Hole	1
8	●	21203	Y-Shape Fracture Plate, 2.0mm, 3 Hole	1
9	●	20606	Y-Shape Fracture Plate, 2.0mm, 4 Hole	1
10	●	20607	Y-Shape Fracture Plate, 2.0mm, 6 Hole	1




Item #	Description	Qty
20478	Cortical Screw, 2.0mm X 8mm	3
20479	Cortical Screw, 2.0mm X 10mm	3
20480	Cortical Screw, 2.0mm X 12mm	3
20481	Cortical Screw, 2.0mm X 14mm	3
20482	Cortical Screw, 2.0mm X 16mm	3
20483	Cortical Screw, 2.0mm X 18mm	3
20484	Cortical Screw, 2.0mm X 20mm	3



Item #	Description	Qty
20471	Cortical Locking Screw, MVA, 2.0mm X 8mm	5
20472	Cortical Locking Screw, MVA, 2.0mm X 10mm	5
20473	Cortical Locking Screw, MVA, 2.0mm X 12mm	5
20474	Cortical Locking Screw, MVA, 2.0mm X 14mm	5
20475	Cortical Locking Screw, MVA, 2.0mm X 16mm	5
20476	Cortical Locking Screw, MVA, 2.0mm X 18mm	5
20477	Cortical Locking Screw, MVA, 2.0mm X 20mm	5

2.0 Screw & Plate Module [continued]



Location	Item #	Description	Qty
11	20636	Locking Drill Guide, 2.0mm	1

Location	Item #	Description	Qty
12	20640	MVA Drill Guide, 2.0mm	1



Location	Item #	Description	Qty
13	20662	Holding & Bending Pin – 2.0	2



Location	Item #	Description	Qty
14	20650	Screwdriver Blade – T6	2






Location	Item #	Description	Qty
15	20647	1.5mm x 95mm Drill Bit (Yellow)	2



Location	Item #	Description	Qty
16	20666	K-Wire With Stop- 1.0mm	2


2.8 Screw Module [ref# 22927]



Item #	Description	Qty
20568	Cortical Screw, 2.8mm X 8mm	3
20569	Cortical Screw, 2.8mm X 10mm	3
20570	Cortical Screw, 2.8mm X 12mm	3
20571	Cortical Screw, 2.8mm X 14mm	3
20572	Cortical Screw, 2.8mm X 16mm	3
20573	Cortical Screw, 2.8mm X 18mm	3
20574	Cortical Screw, 2.8mm X 20mm	3
20575	Cortical Screw, 2.8mm X 22mm	3
20576	Cortical Screw, 2.8mm X 24mm	3
20577	Cortical Screw, 2.8mm X 26mm	3
20578	Cortical Screw, 2.8mm X 28mm	3
20579	Cortical Screw, 2.8mm X 30mm	3

Item #	Description	Qty
20556	Cortical Locking Screw, MVA, 2.8mm X 8mm	5
20557	Cortical Locking Screw, MVA, 2.8mm X 10mm	5
20558	Cortical Locking Screw, MVA, 2.8mm X 12mm	5
20559	Cortical Locking Screw, MVA, 2.8mm X 14mm	5
20560	Cortical Locking Screw, MVA, 2.8mm X 16mm	5
20561	Cortical Locking Screw, MVA, 2.8mm X 18mm	5
20562	Cortical Locking Screw, MVA, 2.8mm X 20mm	5
20563	Cortical Locking Screw, MVA, 2.8mm X 22mm	5
20564	Cortical Locking Screw, MVA, 2.8mm X 24mm	5
20565	Cortical Locking Screw, MVA, 2.8mm X 26mm	5
20566	Cortical Locking Screw, MVA, 2.8mm X 28mm	5
20567	Cortical Locking Screw, MVA, 2.8mm X 30mm	5

2.8 Screw Module [continued]



Location	Item #	Description	Qty
1	20638	Locking Drill Guide, 2.8mm	1



Location	Item #	Description	Qty
2	20641	MVA Drill Guide, 2.5-2.8mm	1



Location	Item #	Description	Qty
3	20663	Holding & Bending Pin – 2.5 – 2.8	2



Location	Item #	Description	Qty
4	20651	Screwdriver Blade – T8	2






Location	Item #	Description	Qty
5	22680	2.0mm x 115mm Drill Bit (Blue)	2



Location	Item #	Description	Qty
6	20940	K-Wire With Stop – 1.2mm	2


3.0/3.5 Screw Module [ref# 22928]




Item #	Description	Qty
20804	Cortical Screw, 3.0mm X 10mm	3
20805	Cortical Screw, 3.0mm X 12mm	3
20806	Cortical Screw, 3.0mm X 14mm	3
20807	Cortical Screw, 3.0mm X 16mm	3
20808	Cortical Screw, 3.0mm X 18mm	3
20809	Cortical Screw, 3.0mm X 20mm	3
20810	Cortical Screw, 3.0mm X 22mm	3
20811	Cortical Screw, 3.0mm X 24mm	3
20812	Cortical Screw, 3.0mm X 26mm	3
20813	Cortical Screw, 3.0mm X 28mm	3
20814	Cortical Screw, 3.0mm X 30mm	3

Item #	Description	Qty
20793	Cortical Locking Screw, MVA, 3.0mm X 10mm	5
20794	Cortical Locking Screw, MVA, 3.0mm X 12mm	5
20795	Cortical Locking Screw, MVA, 3.0mm X 14mm	5
20796	Cortical Locking Screw, MVA, 3.0mm X 16mm	5
20797	Cortical Locking Screw, MVA, 3.0mm X 18mm	5
20798	Cortical Locking Screw, MVA, 3.0mm X 20mm	5
20799	Cortical Locking Screw, MVA, 3.0mm X 22mm	5
20800	Cortical Locking Screw, MVA, 3.0mm X 24mm	5
20801	Cortical Locking Screw, MVA, 3.0mm X 26mm	5
20802	Cortical Locking Screw, MVA, 3.0mm X 28mm	5
20803	Cortical Locking Screw, MVA, 3.0mm X 30mm	5


3.0/3.5 Screw Module [continued]



Item #	Description	Qty
21269	Cortical Screw, Low Head, 3.5mm X 10mm	3
21270	Cortical Screw, Low Head, 3.5mm X 12mm	3
21271	Cortical Screw, Low Head, 3.5mm X 14mm	3
21272	Cortical Screw, Low Head, 3.5mm X 16mm	3
21273	Cortical Screw, Low Head, 3.5mm X 18mm	3
21274	Cortical Screw, Low Head, 3.5mm X 20mm	3
21275	Cortical Screw, Low Head, 3.5mm X 22mm	3
21276	Cortical Screw, Low Head, 3.5mm X 24mm	3
21277	Cortical Screw, Low Head, 3.5mm X 26mm	3
21278	Cortical Screw, Low Head, 3.5mm X 28mm	3
21279	Cortical Screw, Low Head, 3.5mm X 30mm	3
21280	Cortical Screw, Low Head, 3.5mm X 35mm	3
21281	Cortical Screw, Low Head, 3.5mm X 40mm	3
21282	Cortical Screw, Low Head, 3.5mm X 45mm	3
21283	Cortical Screw, Low Head, 3.5mm X 50mm	3
21284	Cortical Screw, Low Head, 3.5mm X 55mm	3



Item #	Description	Qty
20815	Cortical Locking Screw, MVA, 3.5mm X 10mm	5
20816	Cortical Locking Screw, MVA, 3.5mm X 12mm	5
20817	Cortical Locking Screw, MVA, 3.5mm X 14mm	5
20818	Cortical Locking Screw, MVA, 3.5mm X 16mm	5
20819	Cortical Locking Screw, MVA, 3.5mm X 18mm	5
20820	Cortical Locking Screw, MVA, 3.5mm X 20mm	5
20821	Cortical Locking Screw, MVA, 3.5mm X 22mm	5
20822	Cortical Locking Screw, MVA, 3.5mm X 24mm	5
20823	Cortical Locking Screw, MVA, 3.5mm X 26mm	5
20824	Cortical Locking Screw, MVA, 3.5mm X 28mm	5
20825	Cortical Locking Screw, MVA, 3.5mm X 30mm	5
21686	Cortical Locking Screw, MVA, 3.5mm X 32mm	5
21687	Cortical Locking Screw, MVA, 3.5mm X 34mm	5
21688	Cortical Locking Screw, MVA, 3.5mm X 36mm	5
21689	Cortical Locking Screw, MVA, 3.5mm X 38mm	5
20827	Cortical Locking Screw, MVA, 3.5mm X 40mm	5



Location	Item #	Description	Qty
1	20917	Locking Drill Guide, 3.0mm	1
2	20918	Locking Drill Guide, 3.5mm	1



Location	Item #	Description	Qty
3	20921	MVA Drill Guide, 3.0-3.5mm	1

3.0/3.5 Screw Module [continued]



Location	Item #	Description	Qty
4	20739	Holding & Bending Pin – 3.0/3.5	2



Location	Item #	Description	Qty
5	20929	Screwdriver Blade – T10	2



Location	Item #	Description	Qty
6	20924	2.0mm x 115mm Drill Bit (Red)	2
7	20925	2.5mm x 125mm Drill Bit (Orange)	2
8	20926	2.5mm x 150mm Drill Bit (Orange)	2



Location	Item #	Description	Qty
9	20941	K-Wire With Stop – 1.4mm	2


General Instruments




Item #	Description	Qty
20930	Screwdriver Handle – Medium	1
20653	Screwdriver Handle – Small	1



Item #	Description	Qty
20635	Depth Gauge 2.0/2.5/2.8 – 40mm	1
21337	Depth Gauge 3.0/3.5/4.0mm – 50mm	1



Item #	Description	Qty
20643	Double Drill Guide 2.0-2.5mm Screws	1
20922	Double Drill Guide 2.8-3.0mm Screws	1
20923	Double Drill Guide 3.5-4.0mm Screws	1



Item #	Description	Qty
20664	Bending Iron – Left	1
20665	Bending Iron – Right	1



Item #	Description	Qty
20934	Hohman Retractor, 8mm	1
20935	Hohman Retractor, 6mm	1

This content is provided as an educational tool only and is not meant as medical advice in the usage of specific BioPro products. A healthcare professional must use their professional judgment in making any final determinations in product usage and technique. The product's Instructions for Use, should always be reviewed prior to surgery. Postoperative management is patient-specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level or outcomes. This information does not constitute medical, legal, or any other type of professional advice and should not be relied upon as such. It is not to be redistributed, duplicated, or disclosed without the express written consent of BioPro, Inc.



BioPro, Inc.
2929 Lapeer Road, Port Huron, MI 48060, USA
info@bioproimplants.com | 1-810-982-7777
www.bioproimplants.com

