












	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Sterilized using ethylene oxide	Indicates a medical device that has been sterilized using ethylene oxide.
	Sterilized using irradiation	Indicates a medical device that has been sterilized using irradiation.
	Non-sterile	Indicates a medical device that has not been subjected to a sterilization process.
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Serial number	Indicates the manufacturer's serial number so that a specific medical device can be identified.
	Consult instructions for use	Indicates the need for the user to consult the instructions for use.
	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	Date of manufacture	Indicates the date when the medical device was manufactured.
	Manufacturer	Indicates the medical device manufacturer.
	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
	Keep dry	Indicates a medical device that needs protection from moisture.
	Prescription Use Only	Caution: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.