


BIOLOGICALLY ORIENTED PROSTHESES

BIOPRO

Go-EZ™ Headed + HBS® Headless

Cannulated Compression Screw Systems





Unique Cutting Tips

Self-drilling + self-tapping cutting tips for effortless screw advancement.

NEW

HBS Fully Threaded








Fully threaded 3.0mm headless design allows for downsizing without sacrificing thread purchase. The dual thread design (cortical and cancellous) creates compression across the fusion site.



Sterile Packaged

Each screw is individually sterile packed and housed in a protective blister to protect screw threads. This ensures sharp cutting teeth and improves traceability for hospital and patient records.



	Go-EZ Headed Compression Screws							
		Ø	Length	Driver (Hex)	Guide Wire	Drill	Color	Tip*
		2.0mm	6-50mm (2 mm increments)	1.5mm	.035" (0.9mm)	1.7mm	Blue	SDST
		2.5mm	6-50mm (2mm increments)	1.5mm	.035" (0.9mm)	1.7mm	Purple	SDST
		3.0mm	6-40mm (2mm increments)	2.5mm	.045" (1.1mm)	2.0mm	Gold	SDST
		3.5mm	6-50mm (2mm increments)	2.5mm	.045" (1.1mm)	2.0mm	Green	SDST
		4.0mm	6-50mm (2mm increments)	3.0mm	.062" (1.6mm)	2.5mm	Grey	SDST
		4.5mm	6-50mm (2mm increments)	3.0mm	.062" (1.6mm)	2.5mm	Brown	SDST
		6.5mm	30-40mm (2mm Increments) 45-120mm (5mm Increments)	4.0mm	2mm	4.3mm	Rose Gold	SDST

	HBS Headless Compression Screws							
	Ø (Distal)	Ø (Proximal)	Length	Driver (Torx)	Guide Wire	Drill	Color	Tip*
	2.5mm	3.5mm	10-40mm (2mm increments)	T-8	.035" (0.9mm)	1.7mm	Rose Gold	SDST
	3.0mm	3.5mm	10-40mm (2mm increments)	T-8	.035" (0.9mm)	1.7mm	Blue	SDST
	2.5mm	3.5mm	10-30mm (1mm increments)	T-8	.035" (0.9mm)	1.7mm	Purple	ST
	3.0mm	3.9mm	10-30mm (1mm increments)	T-8	.035" (0.9mm)	1.7mm	Grey	ST
	3.0mm	3.9mm	10-30mm (2mm increments)	T-8	.035" (0.9mm)	1.7mm	Silver	SDST

NOTE: Standard surgery kits DO NOT contain all available screw sizes. Full size range is available on special request.

*** SDST = Self-drilling + Self-tapping ST = Self-tapping**

BIOLOGICALLY ORIENTED PROSTHESES

BIOPRO

www.bioproimplants.com

Schedule a surgery today

(810) 982-7777

orders@bioproimplants.com



+1-810-982-7777
info@bioproimplants.com
2929 Lapeer Road
Port Huron, MI 48060



MKT08 02