



Tray Layout Guide

Item # 22675



Top Level

- 1 Plate Module
- 2 Bending Instruments - *Bending Pins, Bending Irons, Hohmann Retractors*
- 3 Plate Module




Bottom Level

- 4 Screw Caddy - *2.0mm MVA Locking and Non-Locking, 2.8mm MVA Locking and Non-Locking, 3.0mm MVA Locking and Non-Locking, 3.5mm MVA Locking, Straight Locking and Non-Locking*
- 5 Drilling Instruments - *Solid Drill Bits, Locking Drill Guides (Tower), MVA Locking Drill Guides (Cone)*
- 6 General Instruments - *Double Drill Guides, Torx Drivers, Olive K-wires, Depth Gauges, Handles*



Note

Torx Drivers are located under double drill guides and Depth Gauges are located under Handles.






| Item # | Description | Qty |
|--------|---|-----|
| 20471 | Cortical Locking Screw, MVA, 2.0mm X 8mm | 3 |
| 20472 | Cortical Locking Screw, MVA, 2.0mm X 10mm | 3 |
| 20473 | Cortical Locking Screw, MVA, 2.0mm X 12mm | 3 |
| 20474 | Cortical Locking Screw, MVA, 2.0mm X 14mm | 3 |
| 20475 | Cortical Locking Screw, MVA, 2.0mm X 16mm | 3 |
| 20476 | Cortical Locking Screw, MVA, 2.0mm X 18mm | 3 |
| 20477 | Cortical Locking Screw, MVA, 2.0mm X 20mm | 3 |

| Item # | Description | Qty |
|--------|------------------------------|-----|
| 20478 | Cortical Screw, 2.0mm X 8mm | 2 |
| 20479 | Cortical Screw, 2.0mm X 10mm | 2 |
| 20480 | Cortical Screw, 2.0mm X 12mm | 2 |
| 20481 | Cortical Screw, 2.0mm X 14mm | 2 |
| 20482 | Cortical Screw, 2.0mm X 16mm | 2 |
| 20483 | Cortical Screw, 2.0mm X 18mm | 2 |
| 20484 | Cortical Screw, 2.0mm X 20mm | 2 |






| Item # | Description | Qty |
|--------|---|-----|
| 20556 | Cortical Locking Screw, MVA, 2.8mm X 8mm | 3 |
| 20557 | Cortical Locking Screw, MVA, 2.8mm X 10mm | 3 |
| 20558 | Cortical Locking Screw, MVA, 2.8mm X 12mm | 3 |
| 20559 | Cortical Locking Screw, MVA, 2.8mm X 14mm | 3 |
| 20560 | Cortical Locking Screw, MVA, 2.8mm X 16mm | 3 |
| 20561 | Cortical Locking Screw, MVA, 2.8mm X 18mm | 3 |
| 20562 | Cortical Locking Screw, MVA, 2.8mm X 20mm | 3 |
| 20563 | Cortical Locking Screw, MVA, 2.8mm X 22mm | 3 |
| 20564 | Cortical Locking Screw, MVA, 2.8mm X 24mm | 3 |
| 20565 | Cortical Locking Screw, MVA, 2.8mm X 26mm | 3 |
| 20566 | Cortical Locking Screw, MVA, 2.8mm X 28mm | 3 |
| 20567 | Cortical Locking Screw, MVA, 2.8mm X 30mm | 3 |

| Item # | Description | Qty |
|--------|------------------------------|-----|
| 20568 | Cortical Screw, 2.8mm X 8mm | 2 |
| 20569 | Cortical Screw, 2.8mm X 10mm | 2 |
| 20570 | Cortical Screw, 2.8mm X 12mm | 2 |
| 20571 | Cortical Screw, 2.8mm X 14mm | 2 |
| 20572 | Cortical Screw, 2.8mm X 16mm | 2 |
| 20573 | Cortical Screw, 2.8mm X 18mm | 2 |
| 20574 | Cortical Screw, 2.8mm X 20mm | 2 |
| 20575 | Cortical Screw, 2.8mm X 22mm | 2 |
| 20576 | Cortical Screw, 2.8mm X 24mm | 2 |
| 20577 | Cortical Screw, 2.8mm X 26mm | 2 |
| 20578 | Cortical Screw, 2.8mm X 28mm | 2 |
| 20579 | Cortical Screw, 2.8mm X 30mm | 2 |




| Item # | Description | Qty |
|--------|---|-----|
| 20793 | Cortical Locking Screw, MVA, 3.0mm X 10mm | 2 |
| 20794 | Cortical Locking Screw, MVA, 3.0mm X 12mm | 2 |
| 20795 | Cortical Locking Screw, MVA, 3.0mm X 14mm | 2 |
| 20796 | Cortical Locking Screw, MVA, 3.0mm X 16mm | 2 |
| 20797 | Cortical Locking Screw, MVA, 3.0mm X 18mm | 2 |
| 20798 | Cortical Locking Screw, MVA, 3.0mm X 20mm | 2 |
| 20799 | Cortical Locking Screw, MVA, 3.0mm X 22mm | 2 |
| 20800 | Cortical Locking Screw, MVA, 3.0mm X 24mm | 2 |
| 20801 | Cortical Locking Screw, MVA, 3.0mm X 26mm | 2 |
| 20802 | Cortical Locking Screw, MVA, 3.0mm X 28mm | 2 |
| 20803 | Cortical Locking Screw, MVA, 3.0mm X 30mm | 2 |

| Item # | Description | Qty |
|--------|------------------------------|-----|
| 20804 | Cortical Screw, 3.0mm X 10mm | 2 |
| 20805 | Cortical Screw, 3.0mm X 12mm | 2 |
| 20806 | Cortical Screw, 3.0mm X 14mm | 2 |
| 20807 | Cortical Screw, 3.0mm X 16mm | 2 |
| 20808 | Cortical Screw, 3.0mm X 18mm | 2 |
| 20809 | Cortical Screw, 3.0mm X 20mm | 2 |
| 20810 | Cortical Screw, 3.0mm X 22mm | 2 |
| 20811 | Cortical Screw, 3.0mm X 24mm | 2 |
| 20812 | Cortical Screw, 3.0mm X 26mm | 2 |
| 20813 | Cortical Screw, 3.0mm X 28mm | 2 |
| 20814 | Cortical Screw, 3.0mm X 30mm | 2 |




| Item # | Description | Qty |
|--------|---|-----|
| 20815 | Cortical Locking Screw, MVA, 3.5mm X 10mm | 3 |
| 20816 | Cortical Locking Screw, MVA, 3.5mm X 12mm | 3 |
| 20817 | Cortical Locking Screw, MVA, 3.5mm X 14mm | 3 |
| 20818 | Cortical Locking Screw, MVA, 3.5mm X 16mm | 3 |
| 20819 | Cortical Locking Screw, MVA, 3.5mm X 18mm | 3 |
| 20820 | Cortical Locking Screw, MVA, 3.5mm X 20mm | 3 |
| 20821 | Cortical Locking Screw, MVA, 3.5mm X 22mm | 3 |
| 20822 | Cortical Locking Screw, MVA, 3.5mm X 24mm | 3 |
| 20823 | Cortical Locking Screw, MVA, 3.5mm X 26mm | 3 |
| 20824 | Cortical Locking Screw, MVA, 3.5mm X 28mm | 3 |
| 20825 | Cortical Locking Screw, MVA, 3.5mm X 30mm | 3 |
| 21686 | Cortical Locking Screw, MVA, 3.5mm X 32mm | 3 |
| 21687 | Cortical Locking Screw, MVA, 3.5mm X 34mm | 3 |
| 21688 | Cortical Locking Screw, MVA, 3.5mm X 36mm | 3 |
| 21689 | Cortical Locking Screw, MVA, 3.5mm X 38mm | 3 |
| 20827 | Cortical Locking Screw, MVA, 3.5mm X 40mm | 3 |

| Item # | Description | Qty |
|--------|-------------------------------------|-----|
| 20852 | Cortical Locking Screw 3.5mm X 10mm | 2 |
| 20853 | Cortical Locking Screw 3.5mm X 12mm | 2 |
| 20854 | Cortical Locking Screw 3.5mm X 14mm | 2 |
| 20855 | Cortical Locking Screw 3.5mm X 16mm | 2 |
| 20856 | Cortical Locking Screw 3.5mm X 18mm | 2 |
| 20857 | Cortical Locking Screw 3.5mm X 20mm | 2 |
| 20858 | Cortical Locking Screw 3.5mm X 22mm | 2 |
| 20859 | Cortical Locking Screw 3.5mm X 24mm | 2 |
| 20860 | Cortical Locking Screw 3.5mm X 26mm | 2 |
| 20861 | Cortical Locking Screw 3.5mm X 28mm | 2 |
| 20862 | Cortical Locking Screw 3.5mm X 30mm | 2 |
| 20863 | Cortical Locking Screw 3.5mm X 35mm | 2 |
| 20864 | Cortical Locking Screw 3.5mm X 40mm | 2 |




| Item # | Description | Qty |
|--------|--|-----|
| 21269 | Cortical Screw, Low Head, 3.5mm X 10mm | 2 |
| 21270 | Cortical Screw, Low Head, 3.5mm X 12mm | 2 |
| 21271 | Cortical Screw, Low Head, 3.5mm X 14mm | 2 |
| 21272 | Cortical Screw, Low Head, 3.5mm X 16mm | 2 |
| 21273 | Cortical Screw, Low Head, 3.5mm X 18mm | 2 |
| 21274 | Cortical Screw, Low Head, 3.5mm X 20mm | 2 |
| 21275 | Cortical Screw, Low Head, 3.5mm X 22mm | 2 |
| 21276 | Cortical Screw, Low Head, 3.5mm X 24mm | 2 |

| Item # | Description | Qty |
|--------|--|-----|
| 21277 | Cortical Screw, Low Head, 3.5mm X 26mm | 2 |
| 21278 | Cortical Screw, Low Head, 3.5mm X 28mm | 2 |
| 21279 | Cortical Screw, Low Head, 3.5mm X 30mm | 2 |
| 21280 | Cortical Screw, Low Head, 3.5mm X 35mm | 2 |
| 21281 | Cortical Screw, Low Head, 3.5mm X 40mm | 2 |
| 21282 | Cortical Screw, Low Head, 3.5mm X 45mm | 2 |
| 21283 | Cortical Screw, Low Head, 3.5mm X 50mm | 2 |
| 21284 | Cortical Screw, Low Head, 3.5mm X 55mm | 2 |




Forefoot Module - 22676

| Location | Screw Color | Item # | Description | Qty |
|----------|-------------|--------|--|-----|
| 1 | •• | 20610 | MTP Fusion Plate, 2.5-3.5mm, Left | 1 |
| 2 | •• | 20611 | MTP Fusion Plate, 2.5-3.5mm, Right | 1 |
| 3 | • | 20628 | Base Open Wedge Plate, 2.5-2.8mm, 0mm Wedge | 2 |
| 4 | • | 20629 | Base Open Wedge Plate, 2.5-2.8mm, 2mm Wedge | 2 |
| 5 | • | 20630 | Base Open Wedge Plate, 2.5-2.8mm, 3mm Wedge | 2 |
| 6 | • | 20631 | Base Open Wedge Plate, 2.5-2.8mm, 4mm Wedge | 2 |
| 7 | • | 20632 | Base Open Wedge Plate, 2.5-2.8mm, 5mm Wedge | 2 |
| 8 | • | 20633 | Base Open Wedge Plate, 2.5-2.8mm, 6mm Wedge | 2 |
| 9 | • | 20605 | Straight Fracture Plate, 2.0mm, 8 Hole | 2 |
| 10 | • | 20604 | Straight Fracture Plate, 2.0mm, 6 Hole | 2 |
| 11 | • | 20602 | Straight Fracture Plate, 2.0mm, 4 Hole | 2 |
| 12 | • | 21207 | Straight Fracture Plate, 2.0mm, 3 Hole | 2 |
| 13 | • | 20607 | Y-Shape Fracture Plate, 2.0mm, 6 Hole | 2 |
| 14 | • | 20606 | Y-Shape Fracture Plate, 2.0mm, 4 Hole | 2 |
| 15 | • | 21203 | Y-Shape Fracture Plate, 2.0mm, 3 Hole | 2 |
| 16 | • | 20609 | T-Shape Fracture Plate, 2.0mm, 6 Hole | 2 |
| 17 | • | 20608 | T-Shape Fracture Plate, 2.0mm, 4 Hole | 2 |
| 18 | • | 21205 | T-Shape Fracture Plate, 2.0mm, 3 Hole | 2 |
| 19 | • | 20621 | L-Shape Fracture Plate, 2.5-2.8mm, 4 Hole, Left | 2 |
| 20 | • | 20623 | L-Shape Fracture Plate, 2.5-2.8mm, 4 Hole, Right | 2 |
| 21 | • | 20620 | L-Shape Fracture Plate, 2.5-2.8mm, 3 Hole, Left | 2 |
| 22 | • | 20622 | L-Shape Fracture Plate, 2.5-2.8mm, 3 Hole, Right | 2 |
| 23 | • | 20626 | Cloverleaf Plate, 2.5-2.8mm, 3 Hole | 2 |
| 24 | • | 20627 | Cloverleaf Plate, 2.5-2.8mm, 4 Hole | 2 |
| 25 | • | 20618 | Y-Shape Fracture Plate, 2.5-2.8mm, 3 Hole | 2 |
| 26 | • | 20634 | Medial TMT-I Fusion Plate, 2.5-2.8mm, 4 Hole | 2 |
| 27 | • | 20619 | Y-Shape Fracture Plate, 2.5-2.8mm, 4 Hole | 2 |
| 28 | • | 20624 | T-Shape Fracture Plate, 2.5-2.8mm, 3 Hole | 2 |
| 29 | • | 20625 | T-Shape Fracture Plate, 2.5-2.8mm, 4 Hole | 2 |
| 30 | • | 20614 | Straight Fracture Plate, 2.5-2.8mm, 4 Hole | 2 |
| 31 | • | 20615 | Straight Fracture Plate, 2.5-2.8mm, 5 Hole | 2 |
| 32 | • | 20616 | Straight Fracture Plate, 2.5-2.8mm, 6 Hole | 2 |
| 33 | • | 20617 | Straight Fracture Plate, 2.5-2.8mm, 7 Hole | 2 |




Midfoot Module - 22671

| Location | Screw Color | Item # | Description | Qty |
|----------|-------------|--------|--|-----|
| 1 | ●● | 20879 | X-Plate, 3.0-3.5mm, Small | 2 |
| 2 | ●● | 20880 | X-Plate, 3.0-3.5mm, Medium | 2 |
| 3 | ●● | 20881 | X-Plate, 3.0-3.5mm, Large | 2 |
| 4 | ●● | 20885 | Rectangular Plate, 3.0-3.5mm, Small | 2 |
| 5 | ●● | 20886 | Rectangular Plate, 3.0-3.5mm, Medium | 2 |
| 6 | ●● | 20887 | Rectangular Plate, 3.0-3.5mm, Large | 2 |
| 7 | ●● | 20874 | L-Shape Fracture Plate, 3.0-3.5mm, 6 Hole, Right | 2 |
| 8 | ●● | 20873 | L-Shape Fracture Plate, 3.0-3.5mm, 4 Hole, Right | 2 |
| 9 | ●● | 21199 | L-Shape Fracture Plate, 3.0-3.5mm, 3 Hole, Right | 2 |
| 10 | ●● | 20878 | T-Shape Fracture Plate, 3.0-3.5mm, 6 Hole | 2 |
| 11 | ●● | 20877 | T-Shape Fracture Plate, 3.0-3.5mm, 4 Hole | 2 |
| 12 | ●● | 21209 | T-Shape Fracture Plate, 3.0-3.5mm, 3 Hole | 2 |
| 13 | ●● | 20870 | Straight Fusion Plate, 3.0-3.5mm, 6 Hole | 2 |
| 14 | ●● | 20869 | Straight Fusion Plate, 3.0-3.5mm, 4 Hole | 2 |
| 15 | ●● | 20868 | Straight Fusion Plate, 3.0-3.5mm, 2 Hole | 2 |
| 16 | ●● | 21201 | Cloverleaf Plate, 3.0-3.5mm, 3 Hole | 2 |
| 17 | ●● | 20875 | Cloverleaf Plate, 3.0-3.5mm, 4 Hole | 2 |
| 18 | ●● | 20876 | Cloverleaf Plate, 3.0-3.5mm, 6 Hole | 2 |
| 19 | ●● | 21197 | L-Shape Fracture Plate, 3.0-3.5mm, 3 Hole, Left | 2 |
| 20 | ●● | 20871 | L-Shape Fracture Plate, 3.0-3.5mm, 4 Hole, Left | 2 |
| 21 | ●● | 20872 | L-Shape Fracture Plate, 3.0-3.5mm, 6 Hole, Left | 2 |




Midfoot Special Module - 22673

| Location | Screw Color | Item # | Description | Qty |
|----------|-------------|--------|---|-----|
| 1 | ●● | 20746 | Dorsal TMT-1 Fusion Plate, 3.0-3.5mm, 0mm Step, Left | 1 |
| 2 | ●● | 20747 | Dorsal TMT-1 Fusion Plate, 3.0-3.5mm, 1mm Step, Left | 1 |
| 3 | ●● | 20748 | Dorsal TMT-1 Fusion Plate, 3.0-3.5mm, 2mm Step, Left | 1 |
| 4 | ●● | 20749 | Dorsal TMT-1 Fusion Plate, 3.0-3.5mm, 3mm Step, Left | 1 |
| 5 | ●● | 20750 | Dorsal TMT-1 Fusion Plate, 3.0-3.5mm, 4mm Step, Left | 1 |
| 6 | ●● | 20751 | Dorsal TMT-1 Fusion Plate, 3.0-3.5mm, 0mm Step, Right | 1 |
| 7 | ●● | 20752 | Dorsal TMT-1 Fusion Plate, 3.0-3.5mm, 1mm Step, Right | 1 |
| 8 | ●● | 20753 | Dorsal TMT-1 Fusion Plate, 3.0-3.5mm, 2mm Step, Right | 1 |
| 9 | ●● | 20754 | Dorsal TMT-1 Fusion Plate, 3.0-3.5mm, 3mm Step, Right | 1 |
| 10 | ●● | 20755 | Dorsal TMT-1 Fusion Plate, 3.0-3.5mm, 4mm Step, Right | 1 |




Hindfoot Module - 22672

| Location | Screw Color | Item # | Description | Qty |
|----------|-------------|--------|--|-----|
| 1 | ●● | 20888 | Evans Plate, 3.0-3.5mm, 4mm Step, Left | 1 |
| 2 | ●● | 20889 | Evans Plate, 3.0-3.5mm, 6mm Step, Left | 1 |
| 3 | ●● | 20890 | Evans Plate, 3.0-3.5mm, 8mm Step, Left | 1 |
| 4 | ●● | 20891 | Evans Plate, 3.0-3.5mm, 10mm Step, Left | 1 |
| 5 | ●● | 20892 | Evans Plate, 3.0-3.5mm, 4mm Step, Right | 1 |
| 6 | ●● | 20893 | Evans Plate, 3.0-3.5mm, 6mm Step, Right | 1 |
| 7 | ●● | 20894 | Evans Plate, 3.0-3.5mm, 8mm Step, Right | 1 |
| 8 | ●● | 20895 | Evans Plate, 3.0-3.5mm, 10mm Step, Right | 1 |
| 9 | ●● | 20896 | Dwyer Plate, 3.0-3.5mm, 4mm Step | 1 |
| 10 | ●● | 20897 | Dwyer Plate, 3.0-3.5mm, 6mm Step | 1 |
| 11 | ●● | 20898 | Dwyer Plate, 3.0-3.5mm, 8mm Step | 1 |
| 12 | ●● | 20899 | Dwyer Plate, 3.0-3.5mm, 10mm Step | 1 |
| 13 | ●● | 20910 | Orif Calcaneal Plate, 3.0-3.5mm, Short | 1 |
| 14 | ●● | 20911 | Orif Calcaneal Plate, 3.0-3.5mm, Medium | 1 |
| 15 | ●● | 20912 | Orif Calcaneal Plate, 3.0-3.5mm, Long | 1 |



Medial Column Module - 22769


| Location | Screw Color | Item # | Description | Qty |
|----------|-------------|--------|--|-----|
| 1 | ●● | 20909 | Medial Distal Column Plate, 3.0-3.5mm, Short | 1 |
| 2 | ●● | 20908 | Medial Distal Column Plate, 3.0-3.5mm, Long | 1 |
| 3 | ●● | 20905 | Medial 3 Column Plate, 3.0-3.5mm, Short, Left | 1 |
| 4 | ●● | 20904 | Medial 3 Column Plate, 3.0-3.5mm, Short, Right | 1 |
| 5 | ●● | 20907 | Medial 3 Column Plate, 3.0-3.5mm, Long, Left | 1 |
| 6 | ●● | 20906 | Medial 3 Column Plate, 3.0-3.5mm, Long, Right | 1 |
| 7 | ●● | 20901 | Medial 4 Column Plate, 3.0-3.5mm, Short, Left | 1 |
| 8 | ●● | 20900 | Medial 4 Column Plate, 3.0-3.5mm, Short, Right | 1 |
| 9 | ●● | 20903 | Medial 4 Column Plate, 3.0-3.5mm, Long, Left | 1 |
| 10 | ●● | 20902 | Medial 4 Column Plate, 3.0-3.5mm, Long, Right | 1 |



| Location | Screw Color | Item # | Description | Qty |
|----------|-------------|--------|-------------------------------------|-----|
| 1 | ● | 20650 | Screwdriver Blade – T6 | 2 |
| 2 | ● | 20651 | Screwdriver Blade – T8 | 2 |
| 3 | ●● | 20929 | Screwdriver Blade – T10 | 2 |
| 4 | ● | 20643 | Double Drill Guide 2.0-2.5mm Screws | 1 |
| 5 | ●● | 20922 | Double Drill Guide 2.8-3.0mm Screws | 1 |
| 6 | ● | 20923 | Double Drill Guide 3.5-4.0mm Screws | 1 |
| 7 | ●● | 20635 | Depth Gauge 2.0/2.5/2.8 – 40mm | 1 |
| 8 | ●● | 21337 | Depth Gauge 3.0/3.5/4.0mm – 50mm | 1 |
| 9 | N/A | 20930 | Screwdriver Handle – Medium | 1 |
| 10 | N/A | 20653 | Screwdriver Handle – Small | 1 |
| 11 | ● | 20666 | K-Wire With Stop- 1.0mm | 2 |
| 12 | ● | 20940 | K-Wire With Stop – 1.2mm | 2 |
| 13 | ●● | 20941 | K-Wire With Stop – 1.4mm | 2 |

Caution


Always ensure that the color marking on the instrument matches the color of the screw.



| Location | Screw Color | Item # | Description | Qty |
|----------|-------------|--------|----------------------------|-----|
| 1 | ● | 20647 | 1.5mm x 95mm Drill Bit | 2 |
| 2 | ● | 22680 | 2.0mm x 115mm Drill Bit | 2 |
| 3 | ● | 20924 | 2.0mm x 115mm Drill Bit | 2 |
| 4 | ● | 20925 | 2.5mm x 125mm Drill Bit | 2 |
| 5 | ● | 20926 | 2.5mm x 150mm Drill Bit | 2 |
| 6 | ● | 20640 | MVA Drill Guide, 2.0mm | 1 |
| 7 | ● | 20641 | MVA Drill Guide, 2.5-2.8mm | 1 |
| 8 | ●● | 20921 | MVA Drill Guide, 3.0-3.5mm | 1 |
| 9 | ● | 20636 | Locking Drill Guide, 2.0mm | 1 |
| 10 | ● | 20638 | Locking Drill Guide, 2.8mm | 1 |
| 11 | ● | 20917 | Locking Drill Guide, 3.0mm | 1 |
| 12 | ● | 20918 | Locking Drill Guide, 3.5mm | 1 |

Caution

Always ensure that the color marking on the instrument matches the color of the screw.



| Location | Screw Color | Item # | Description | Qty |
|----------|-------------|--------|---------------------------------------|-----|
| 1 | ● | 20662 | Holding & Bending Pin – 2.0mm | 2 |
| 2 | ● | 20663 | Holding & Bending Pin – 2.5 – 2.8mm | 2 |
| 3 | ●● | 20739 | Holding & Bending Pin – 3.0/3.5/4.0mm | 2 |
| 4 | N/A | 20664 | Bending Iron – Left | 1 |
| 5 | N/A | 20665 | Bending Iron – Right | 1 |
| 6 | N/A | 20934 | Hohman Retractor, 8mm | 1 |
| 7 | N/A | 20935 | Hohman Retractor, 6mm | 1 |

Caution

Always ensure that the color marking on the instrument matches the color of the screw.

To place an order email orders@bioproimplants.com or call 810-982-7777



+1-810-982-7777
info@bioproimplants.com
 2929 Lapeer Road
 Port Huron, MI 48060

