



## Tray Layout Guide

The BioPro Foot Plating System tray utilizes a modular design allowing the user to customize the tray, thereby reducing the footprint in the operating room.

- Base Instrument Tray [ref #22925]
- 3.0/3.5 Plate Module [ref #22931]
- 2.8 Plate Module [ref #22930]
- Medial Column Plate Module [ref #22932]
- Calcaneal Plate Module [ref #22933]
- MTP Plate Module [ref #22929]
- 2.0 Screw + Plate Module [ref #22926]
- 2.8 Screw Module [ref #22927]
- 3.0/3.5 Screw Module [ref #22928]
- Osteotomy Plate Module [ref #22934]




The Foot Plating System Main Tray [ref# 22925] is constructed with 3 levels. Each module is to be ordered separately.




#### Top Level

- 1 Plate Modules (ordered separately)
  - Medial Column Plate Module [ref #22932]
  - Calcaneal Plate Module [ref #22933]
- 2 Plate Module (ordered separately)
  - 3.0/3.5 Plate Module [ref #22931]
  - 2.8 Plate Module [ref #22930]
- 3 Drivers
- 4 Depth Gauges



#### Middle Level


- 5 3.0/3.5 Screw Module [ref #22928] (ordered separately)
- 6 2.8 Screw Module [ref #22927] (ordered separately)
- 7 2.0 Screw & Plate Module [ref #22926] (ordered separately)
- 8 MTP Module [ref #22929] (ordered separately)



#### Bottom Level

- 9 Osteotomy Plate Module [ref #22934] (ordered separately)
- 10 Bending Irons
- 11 Double Drill Guides
- 12 Hohman Retractors
- 13 Miscellaneous Instruments (ordered separately)


## 2.8 Plate Module [ref# 22930]



**2.8 Plate Module - 22930**

| Location | Screw Color | Item # | Description                                      | Qty |
|----------|-------------|--------|--------------------------------------------------|-----|
| 1        | •           | 20614  | Straight Fracture Plate, 2.5-2.8mm, 4 Hole       | 1   |
| 2        | •           | 20615  | Straight Fracture Plate, 2.5-2.8mm, 5 Hole       | 1   |
| 3        | •           | 20616  | Straight Fracture Plate, 2.5-2.8mm, 6 Hole       | 1   |
| 4        | •           | 20617  | Straight Fracture Plate, 2.5-2.8mm, 7 Hole       | 1   |
| 5        | •           | 20624  | T-Shape Fracture Plate, 2.5-2.8mm, 3 Hole        | 1   |
| 6        | •           | 20625  | T-Shape Fracture Plate, 2.5-2.8mm, 4 Hole        | 1   |
| 7        | •           | 20622  | L-Shape Fracture Plate, 2.5-2.8mm, 3 Hole, Right | 1   |
| 8        | •           | 20623  | L-Shape Fracture Plate, 2.5-2.8mm, 4 Hole, Right | 1   |
| 9        | •           | 20620  | L-Shape Fracture Plate, 2.5-2.8mm, 3 Hole, Left  | 1   |
| 10       | •           | 20621  | L-Shape Fracture Plate, 2.5-2.8mm, 4 Hole, Left  | 1   |
| 11       | •           | 20756  | Medial TMT-1 Fusion Plate, 2.8mm, 4 Hole         | 1   |
| 12       | •           | 20618  | Y-Shape Fracture Plate, 2.5-2.8mm, 3 Hole        | 1   |
| 13       | •           | 20619  | Y-Shape Fracture Plate, 2.5-2.8mm, 4 Hole        | 1   |
| 14       | •           | 20626  | Cloverleaf Plate, 2.5-2.8mm, 3 Hole              | 1   |
| 15       | •           | 20627  | Cloverleaf Plate, 2.5-2.8mm, 4 Hole              | 1   |
| 16       | •           | 20628  | Open Wedge Plate, 2.5-2.8mm, 0mm Wedge           | 1   |
| 17       | •           | 20629  | Open Wedge Plate, 2.5-2.8mm, 2mm Wedge           | 1   |
| 18       | •           | 20630  | Open Wedge Plate, 2.5-2.8mm, 3mm Wedge           | 1   |
| 19       | •           | 20631  | Open Wedge Plate, 2.5-2.8mm, 4mm Wedge           | 1   |
| 20       | •           | 20632  | Open Wedge Plate, 2.5-2.8mm, 5mm Wedge           | 1   |
| 21       | •           | 20633  | Open Wedge Plate, 2.5-2.8mm, 6mm Wedge           | 1   |


## 3.0-3.5 Plate Module [ref# 22931]



**3.0-3.5 Plate Module - 22931**

| Location | Screw Color | Item # | Description                                      | Qty |
|----------|-------------|--------|--------------------------------------------------|-----|
| 1        | ● ●         | 20868  | Straight Fusion Plate, 3.0-3.5mm, 2 Hole         | 1   |
| 2        | ● ●         | 20869  | Straight Fusion Plate, 3.0-3.5mm, 4 Hole         | 1   |
| 3        | ● ●         | 20870  | Straight Fusion Plate, 3.0-3.5mm, 6 Hole         | 1   |
| 4        | ● ●         | 21209  | T-Shape Fracture Plate, 3.0-3.5mm, 3 Hole        | 1   |
| 5        | ● ●         | 20877  | T-Shape Fracture Plate, 3.0-3.5mm, 4 Hole        | 1   |
| 6        | ● ●         | 20878  | T-Shape Fracture Plate, 3.0-3.5mm, 6 Hole        | 1   |
| 7        | ● ●         | 21197  | L-Shape Fracture Plate, 3.0-3.5mm, 3 Hole, Left  | 1   |
| 8        | ● ●         | 20871  | L-Shape Fracture Plate, 3.0-3.5mm, 4 Hole, Left  | 1   |
| 9        | ● ●         | 20872  | L-Shape Fracture Plate, 3.0-3.5mm, 6 Hole, Left  | 1   |
| 10       | ● ●         | 21199  | L-Shape Fracture Plate, 3.0-3.5mm, 3 Hole, Right | 1   |
| 11       | ● ●         | 20873  | L-Shape Fracture Plate, 3.0-3.5mm, 4 Hole, Right | 1   |
| 12       | ● ●         | 20874  | L-Shape Fracture Plate, 3.0-3.5mm, 6 Hole, Right | 1   |
| 13       | ● ●         | 21201  | Cloverleaf Plate, 3.0-3.5mm, 3 Hole              | 1   |
| 14       | ● ●         | 20875  | Cloverleaf Plate, 3.0-3.5mm, 4 Hole              | 1   |
| 15       | ● ●         | 20876  | Cloverleaf Plate, 3.0-3.5mm, 6 Hole              | 1   |
| 16       | ● ●         | 20879  | X-Plate, 3.0-4.0mm, Small                        | 1   |
| 17       | ● ●         | 20880  | X-Plate, 3.0-4.0mm, Medium                       | 1   |
| 18       | ● ●         | 20881  | X-Plate, 3.0-4.0mm, Large                        | 1   |
| 19       | ● ●         | 20885  | Rectangular Plate, 3.0-3.5mm, Small              | 1   |
| 20       | ● ●         | 20886  | Rectangular Plate, 3.0-3.5mm, Medium             | 1   |
| 21       | ● ●         | 20887  | Rectangular Plate, 3.0-3.5mm, Large              | 1   |


MTP Plate Module [ref# 22929]



| MTP Plate Module - 22929 |             |        |                                        |     |
|--------------------------|-------------|--------|----------------------------------------|-----|
| Location                 | Screw Color | Item # | Description                            | Qty |
| 1                        | ● ● ●       | 20611  | MTP Fusion Plate, 2.5-3.5mm, Right     | 1   |
| 2                        | ● ● ●       | 20610  | MTP Fusion Plate, 2.5-3.5mm, Left      | 1   |
| 3                        | ● ● ●       | 20613  | MTP Fusion Rev Plate, 2.8-3.5mm, Right | 1   |
| 4                        | ● ● ●       | 20612  | MTP Fusion Rev Plate, 2.8-3.5mm, Left  | 1   |
| 5                        |             | 20671  | MTP Fusion Reamer, Convex, 16mm        | 1   |
| 6                        |             | 20672  | MTP Fusion Reamer, Convex, 18mm        | 1   |
| 7                        |             | 20673  | MTP Fusion Reamer, Convex, 21mm        | 1   |
| 8                        |             | 20674  | MTP Fusion Reamer, Convex, 24mm        | 1   |
| 9                        |             | 20677  | MTP Fusion Reamer, Concave, 16mm       | 1   |
| 10                       |             | 20678  | MTP Fusion Reamer, Concave, 18mm       | 1   |
| 11                       |             | 20679  | MTP Fusion Reamer, Concave, 21mm       | 1   |
| 12                       |             | 20680  | MTP Fusion Reamer, Concave, 24mm       | 1   |
| 13                       |             | 17232  | K-Wire .062" x 4" Non-Sterile          | 4   |




# Osteotomy Plate Module [ref# 22934]



**Osteotomy Module - 22934**


| Location | Screw Color | Item # | Description                              | Qty |
|----------|-------------|--------|------------------------------------------|-----|
| 1        | ● ●         | 20888  | Evans Plate, 3.0-3.5mm, 4mm Step, Left   | 1   |
| 2        | ● ●         | 20889  | Evans Plate, 3.0-3.5mm, 6mm Step, Left   | 1   |
| 3        | ● ●         | 20890  | Evans Plate, 3.0-3.5mm, 8mm Step, Left   | 1   |
| 4        | ● ●         | 20891  | Evans Plate, 3.0-3.5mm, 10mm Step, Left  | 1   |
| 5        | ● ●         | 20892  | Evans Plate, 3.0-3.5mm, 4mm Step, Right  | 1   |
| 6        | ● ●         | 20893  | Evans Plate, 3.0-3.5mm, 6mm Step, Right  | 1   |
| 7        | ● ●         | 20894  | Evans Plate, 3.0-3.5mm, 8mm Step, Right  | 1   |
| 8        | ● ●         | 20895  | Evans Plate, 3.0-3.5mm, 10mm Step, Right | 1   |
| 9        | ● ●         | 20896  | Dwyer Plate, 3.0-3.5mm, 4mm Step         | 1   |
| 10       | ● ●         | 20897  | Dwyer Plate, 3.0-3.5mm, 6mm Step         | 1   |
| 11       | ● ●         | 20898  | Dwyer Plate, 3.0-3.5mm, 8mm Step         | 1   |
| 12       | ● ●         | 20899  | Dwyer Plate, 3.0-3.5mm, 10mm Step        | 1   |

Medial Column Plate Module [ref# 22932]



| Medial Column Module - 22932 |             |        |                                                |     |
|------------------------------|-------------|--------|------------------------------------------------|-----|
| Location                     | Screw Color | Item # | Description                                    | Qty |
| 1                            | ● ●         | 20909  | Medial Distal Column Plate, 3.0-3.5mm, Short   | 1   |
| 2                            | ● ●         | 20908  | Medial Distal Column Plate, 3.0-3.5mm, Long    | 1   |
| 3                            | ● ●         | 20905  | Medial 3 Column Plate, 3.0-3.5mm, Short, Left  | 1   |
| 4                            | ● ●         | 20904  | Medial 3 Column Plate, 3.0-3.5mm, Short, Right | 1   |
| 5                            | ● ●         | 20907  | Medial 3 Column Plate, 3.0-3.5mm, Long, Left   | 1   |
| 6                            | ● ●         | 20906  | Medial 3 Column Plate, 3.0-3.5mm, Long, Right  | 1   |
| 7                            | ● ●         | 20901  | Medial 4 Column Plate, 3.0-3.5mm, Short, Left  | 1   |
| 8                            | ● ●         | 20900  | Medial 4 Column Plate, 3.0-3.5mm, Short, Right | 1   |
| 9                            | ● ●         | 20903  | Medial 4 Column Plate, 3.0-3.5mm, Long, Left   | 1   |
| 10                           | ● ●         | 20902  | Medial 4 Column Plate, 3.0-3.5mm, Long, Right  | 1   |


Calcaneal Plate Module [ref# 22933]



| Calcaneal Plate Module - 22933 |             |        |                                         |     |
|--------------------------------|-------------|--------|-----------------------------------------|-----|
| Location                       | Screw Color | Item # | Description                             | Qty |
| 1                              | ● ●         | 20910  | ORIF Calcaneal Plate, 3.0-3.5mm, Short  | 1   |
| 2                              | ● ●         | 20911  | ORIF Calcaneal Plate, 3.0-3.5mm, Medium | 1   |
| 3                              | ● ●         | 20912  | ORIF Calcaneal Plate, 3.0-3.5mm, Long   | 1   |



2.0 Screw & Plate Module [ref# 22926]



2.0 Screw-Plate Module - 22926

| Location | Screw Color | Item # | Description                            | Qty |
|----------|-------------|--------|----------------------------------------|-----|
| 1        | ●           | 21207  | Straight Fracture Plate, 2.0mm, 3 Hole | 1   |
| 2        | ●           | 20602  | Straight Fracture Plate, 2.0mm, 4 Hole | 1   |
| 3        | ●           | 20604  | Straight Fracture Plate, 2.0mm, 6 Hole | 1   |
| 4        | ●           | 20605  | Straight Fracture Plate, 2.0mm, 8 Hole | 1   |
| 5        | ●           | 21205  | T-Shape Fracture Plate, 2.0mm, 3 Hole  | 1   |
| 6        | ●           | 20608  | T-Shape Fracture Plate, 2.0mm, 4 Hole  | 1   |
| 7        | ●           | 20609  | T-Shape Fracture Plate, 2.0mm, 6 Hole  | 1   |
| 8        | ●           | 21203  | Y-Shape Fracture Plate, 2.0mm, 3 Hole  | 1   |
| 9        | ●           | 20606  | Y-Shape Fracture Plate, 2.0mm, 4 Hole  | 1   |
| 10       | ●           | 20607  | Y-Shape Fracture Plate, 2.0mm, 6 Hole  | 1   |




| Item # | Description                  | Qty |
|--------|------------------------------|-----|
| 20478  | Cortical Screw, 2.0mm X 8mm  | 3   |
| 20479  | Cortical Screw, 2.0mm X 10mm | 3   |
| 20480  | Cortical Screw, 2.0mm X 12mm | 3   |
| 20481  | Cortical Screw, 2.0mm X 14mm | 3   |
| 20482  | Cortical Screw, 2.0mm X 16mm | 3   |
| 20483  | Cortical Screw, 2.0mm X 18mm | 3   |
| 20484  | Cortical Screw, 2.0mm X 20mm | 3   |



| Item # | Description                               | Qty |
|--------|-------------------------------------------|-----|
| 20471  | Cortical Locking Screw, MVA, 2.0mm X 8mm  | 5   |
| 20472  | Cortical Locking Screw, MVA, 2.0mm X 10mm | 5   |
| 20473  | Cortical Locking Screw, MVA, 2.0mm X 12mm | 5   |
| 20474  | Cortical Locking Screw, MVA, 2.0mm X 14mm | 5   |
| 20475  | Cortical Locking Screw, MVA, 2.0mm X 16mm | 5   |
| 20476  | Cortical Locking Screw, MVA, 2.0mm X 18mm | 5   |
| 20477  | Cortical Locking Screw, MVA, 2.0mm X 20mm | 5   |

2.0 Screw & Plate Module [continued]



| Location | Item # | Description                | Qty |
|----------|--------|----------------------------|-----|
| 11       | 20636  | Locking Drill Guide, 2.0mm | 1   |

| Location | Item # | Description            | Qty |
|----------|--------|------------------------|-----|
| 12       | 20640  | MVA Drill Guide, 2.0mm | 1   |



| Location | Item # | Description                 | Qty |
|----------|--------|-----------------------------|-----|
| 13       | 20662  | Holding & Bending Pin – 2.0 | 2   |



| Location | Item # | Description            | Qty |
|----------|--------|------------------------|-----|
| 14       | 20650  | Screwdriver Blade – T6 | 2   |






| Location | Item # | Description                     | Qty |
|----------|--------|---------------------------------|-----|
| 15       | 20647  | 1.5mm x 95mm Drill Bit (Yellow) | 2   |



| Location | Item # | Description             | Qty |
|----------|--------|-------------------------|-----|
| 16       | 20666  | K-Wire With Stop- 1.0mm | 2   |

2.8 Screw Module [ref# 22927]



| Item # | Description                  | Qty |
|--------|------------------------------|-----|
| 20568  | Cortical Screw, 2.8mm X 8mm  | 3   |
| 20569  | Cortical Screw, 2.8mm X 10mm | 3   |
| 20570  | Cortical Screw, 2.8mm X 12mm | 3   |
| 20571  | Cortical Screw, 2.8mm X 14mm | 3   |
| 20572  | Cortical Screw, 2.8mm X 16mm | 3   |
| 20573  | Cortical Screw, 2.8mm X 18mm | 3   |
| 20574  | Cortical Screw, 2.8mm X 20mm | 3   |
| 20575  | Cortical Screw, 2.8mm X 22mm | 3   |
| 20576  | Cortical Screw, 2.8mm X 24mm | 3   |
| 20577  | Cortical Screw, 2.8mm X 26mm | 3   |
| 20578  | Cortical Screw, 2.8mm X 28mm | 3   |
| 20579  | Cortical Screw, 2.8mm X 30mm | 3   |

| Item # | Description                               | Qty |
|--------|-------------------------------------------|-----|
| 20556  | Cortical Locking Screw, MVA, 2.8mm X 8mm  | 5   |
| 20557  | Cortical Locking Screw, MVA, 2.8mm X 10mm | 5   |
| 20558  | Cortical Locking Screw, MVA, 2.8mm X 12mm | 5   |
| 20559  | Cortical Locking Screw, MVA, 2.8mm X 14mm | 5   |
| 20560  | Cortical Locking Screw, MVA, 2.8mm X 16mm | 5   |
| 20561  | Cortical Locking Screw, MVA, 2.8mm X 18mm | 5   |
| 20562  | Cortical Locking Screw, MVA, 2.8mm X 20mm | 5   |
| 20563  | Cortical Locking Screw, MVA, 2.8mm X 22mm | 5   |
| 20564  | Cortical Locking Screw, MVA, 2.8mm X 24mm | 5   |
| 20565  | Cortical Locking Screw, MVA, 2.8mm X 26mm | 5   |
| 20566  | Cortical Locking Screw, MVA, 2.8mm X 28mm | 5   |
| 20567  | Cortical Locking Screw, MVA, 2.8mm X 30mm | 5   |

## 2.8 Screw Module [continued]



| Location | Item # | Description                | Qty |
|----------|--------|----------------------------|-----|
| 1        | 20638  | Locking Drill Guide, 2.8mm | 1   |



| Location | Item # | Description                | Qty |
|----------|--------|----------------------------|-----|
| 2        | 20641  | MVA Drill Guide, 2.5-2.8mm | 1   |



| Location | Item # | Description                       | Qty |
|----------|--------|-----------------------------------|-----|
| 3        | 20663  | Holding & Bending Pin – 2.5 – 2.8 | 2   |



| Location | Item # | Description            | Qty |
|----------|--------|------------------------|-----|
| 4        | 20651  | Screwdriver Blade – T8 | 2   |






| Location | Item # | Description                    | Qty |
|----------|--------|--------------------------------|-----|
| 5        | 22680  | 2.0mm x 115mm Drill Bit (Blue) | 2   |



| Location | Item # | Description              | Qty |
|----------|--------|--------------------------|-----|
| 6        | 20940  | K-Wire With Stop – 1.2mm | 2   |


3.0/3.5 Screw Module [ref# 22928]




| Item # | Description                  | Qty |
|--------|------------------------------|-----|
| 20804  | Cortical Screw, 3.0mm X 10mm | 3   |
| 20805  | Cortical Screw, 3.0mm X 12mm | 3   |
| 20806  | Cortical Screw, 3.0mm X 14mm | 3   |
| 20807  | Cortical Screw, 3.0mm X 16mm | 3   |
| 20808  | Cortical Screw, 3.0mm X 18mm | 3   |
| 20809  | Cortical Screw, 3.0mm X 20mm | 3   |
| 20810  | Cortical Screw, 3.0mm X 22mm | 3   |
| 20811  | Cortical Screw, 3.0mm X 24mm | 3   |
| 20812  | Cortical Screw, 3.0mm X 26mm | 3   |
| 20813  | Cortical Screw, 3.0mm X 28mm | 3   |
| 20814  | Cortical Screw, 3.0mm X 30mm | 3   |

| Item # | Description                               | Qty |
|--------|-------------------------------------------|-----|
| 20793  | Cortical Locking Screw, MVA, 3.0mm X 10mm | 5   |
| 20794  | Cortical Locking Screw, MVA, 3.0mm X 12mm | 5   |
| 20795  | Cortical Locking Screw, MVA, 3.0mm X 14mm | 5   |
| 20796  | Cortical Locking Screw, MVA, 3.0mm X 16mm | 5   |
| 20797  | Cortical Locking Screw, MVA, 3.0mm X 18mm | 5   |
| 20798  | Cortical Locking Screw, MVA, 3.0mm X 20mm | 5   |
| 20799  | Cortical Locking Screw, MVA, 3.0mm X 22mm | 5   |
| 20800  | Cortical Locking Screw, MVA, 3.0mm X 24mm | 5   |
| 20801  | Cortical Locking Screw, MVA, 3.0mm X 26mm | 5   |
| 20802  | Cortical Locking Screw, MVA, 3.0mm X 28mm | 5   |
| 20803  | Cortical Locking Screw, MVA, 3.0mm X 30mm | 5   |


3.0/3.5 Screw Module [continued]



| Item # | Description                            | Qty |
|--------|----------------------------------------|-----|
| 21269  | Cortical Screw, Low Head, 3.5mm X 10mm | 3   |
| 21270  | Cortical Screw, Low Head, 3.5mm X 12mm | 3   |
| 21271  | Cortical Screw, Low Head, 3.5mm X 14mm | 3   |
| 21272  | Cortical Screw, Low Head, 3.5mm X 16mm | 3   |
| 21273  | Cortical Screw, Low Head, 3.5mm X 18mm | 3   |
| 21274  | Cortical Screw, Low Head, 3.5mm X 20mm | 3   |
| 21275  | Cortical Screw, Low Head, 3.5mm X 22mm | 3   |
| 21276  | Cortical Screw, Low Head, 3.5mm X 24mm | 3   |
| 21277  | Cortical Screw, Low Head, 3.5mm X 26mm | 3   |
| 21278  | Cortical Screw, Low Head, 3.5mm X 28mm | 3   |
| 21279  | Cortical Screw, Low Head, 3.5mm X 30mm | 3   |
| 21280  | Cortical Screw, Low Head, 3.5mm X 35mm | 3   |
| 21281  | Cortical Screw, Low Head, 3.5mm X 40mm | 3   |
| 21282  | Cortical Screw, Low Head, 3.5mm X 45mm | 3   |
| 21283  | Cortical Screw, Low Head, 3.5mm X 50mm | 3   |
| 21284  | Cortical Screw, Low Head, 3.5mm X 55mm | 3   |



| Item # | Description                               | Qty |
|--------|-------------------------------------------|-----|
| 20815  | Cortical Locking Screw, MVA, 3.5mm X 10mm | 5   |
| 20816  | Cortical Locking Screw, MVA, 3.5mm X 12mm | 5   |
| 20817  | Cortical Locking Screw, MVA, 3.5mm X 14mm | 5   |
| 20818  | Cortical Locking Screw, MVA, 3.5mm X 16mm | 5   |
| 20819  | Cortical Locking Screw, MVA, 3.5mm X 18mm | 5   |
| 20820  | Cortical Locking Screw, MVA, 3.5mm X 20mm | 5   |
| 20821  | Cortical Locking Screw, MVA, 3.5mm X 22mm | 5   |
| 20822  | Cortical Locking Screw, MVA, 3.5mm X 24mm | 5   |
| 20823  | Cortical Locking Screw, MVA, 3.5mm X 26mm | 5   |
| 20824  | Cortical Locking Screw, MVA, 3.5mm X 28mm | 5   |
| 20825  | Cortical Locking Screw, MVA, 3.5mm X 30mm | 5   |
| 21686  | Cortical Locking Screw, MVA, 3.5mm X 32mm | 5   |
| 21687  | Cortical Locking Screw, MVA, 3.5mm X 34mm | 5   |
| 21688  | Cortical Locking Screw, MVA, 3.5mm X 36mm | 5   |
| 21689  | Cortical Locking Screw, MVA, 3.5mm X 38mm | 5   |
| 20827  | Cortical Locking Screw, MVA, 3.5mm X 40mm | 5   |



| Location | Item # | Description                | Qty |
|----------|--------|----------------------------|-----|
| 1        | 20917  | Locking Drill Guide, 3.0mm | 1   |
| 2        | 20918  | Locking Drill Guide, 3.5mm | 1   |



| Location | Item # | Description                | Qty |
|----------|--------|----------------------------|-----|
| 3        | 20921  | MVA Drill Guide, 3.0-3.5mm | 1   |



3.0/3.5 Screw Module [continued]



| Location | Item # | Description                     | Qty |
|----------|--------|---------------------------------|-----|
| 4        | 20739  | Holding & Bending Pin – 3.0/3.5 | 2   |



| Location | Item # | Description             | Qty |
|----------|--------|-------------------------|-----|
| 5        | 20929  | Screwdriver Blade – T10 | 2   |



| Location | Item # | Description                      | Qty |
|----------|--------|----------------------------------|-----|
| 6        | 20924  | 2.0mm x 115mm Drill Bit (Red)    | 2   |
| 7        | 20925  | 2.5mm x 125mm Drill Bit (Orange) | 2   |
| 8        | 20926  | 2.5mm x 150mm Drill Bit (Orange) | 2   |



| Location | Item # | Description              | Qty |
|----------|--------|--------------------------|-----|
| 9        | 20941  | K-Wire With Stop – 1.4mm | 2   |


# General Instruments




| Item # | Description                 | Qty |
|--------|-----------------------------|-----|
| 20930  | Screwdriver Handle – Medium | 1   |
| 20653  | Screwdriver Handle – Small  | 1   |



| Item # | Description                      | Qty |
|--------|----------------------------------|-----|
| 20635  | Depth Gauge 2.0/2.5/2.8 – 40mm   | 1   |
| 21337  | Depth Gauge 3.0/3.5/4.0mm – 50mm | 1   |



| Item # | Description                         | Qty |
|--------|-------------------------------------|-----|
| 20643  | Double Drill Guide 2.0-2.5mm Screws | 1   |
| 20922  | Double Drill Guide 2.8-3.0mm Screws | 1   |
| 20923  | Double Drill Guide 3.5-4.0mm Screws | 1   |



| Item # | Description          | Qty |
|--------|----------------------|-----|
| 20664  | Bending Iron – Left  | 1   |
| 20665  | Bending Iron – Right | 1   |



| Item # | Description           | Qty |
|--------|-----------------------|-----|
| 20934  | Hohman Retractor, 8mm | 1   |
| 20935  | Hohman Retractor, 6mm | 1   |

This content is provided as an educational tool only and is not meant as medical advice in the usage of specific BioPro products. A healthcare professional must use their professional judgment in making any final determinations in product usage and technique. The product's Instructions for Use, should always be reviewed prior to surgery. Postoperative management is patient-specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level or outcomes. This information does not constitute medical, legal, or any other type of professional advice and should not be relied upon as such. It is not to be redistributed, duplicated, or disclosed without the express written consent of BioPro, Inc.