


Foot & Ankle

Our comprehensive Foot and Ankle portfolio offers straightforward solutions to address complex conditions including Hallux Rigidus, Hallux Valgus, Flatfoot, Hammertoes, Fractures, and more.





Accu-Cut® Osteotomy Guide System

Metatarsal osteotomy guides in an all-inclusive sterile package provide the surgeon with everything needed to perform the desired osteotomy.



MPJ Hemi Implants

A proven fusion alternative backed by 70+ years of clinical data. The procedure offers restoration of motion and pain relief, with a minimal bone resection technique.




Forefoot-Hindfoot Plating System

A comprehensive line of multi-variable angle locking plates to address a variety of procedures. The system features 2.0, 2.8, 3.0 and 3.5mm screws.



Digital Compression Screw

A 1.5 and 1.8mm solid, stainless steel screw specifically designed to address digital fusions.



Tendon Anchor System™

A toothed titanium anchor and screw construct designed for fast and secure tendon attachment.



Horizon® Subtalar

A threaded arthroereisis implant with a smooth lateral surface to address flexible flat foot with a quick and reversible procedure.



Cannulated Compression Screws

A comprehensive line of titanium cannulated screws including headed and headless screws. Go-EZ® headed screws range from 2.0 to 6.5mm while headless designs are available in 2.5 or 3.0mm.



Memory Staple

A 2-leg nitinol memory-alloy staple available in 17 standard sizes. The unique s-bridge staple activates at patient body temperature to provide dynamic compression.

This content is provided as an educational tool only and is not meant as medical advice in the usage of specific BioPro products. A healthcare professional must use their professional judgment in making any final determinations in product usage and technique. The product's Instructions for Use, should always be reviewed prior to surgery. Postoperative management is patient-specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level or outcomes. This information does not constitute medical, legal, or any other type of professional advice and should not be relied upon as such. It is not to be redistributed, duplicated, or disclosed without the express written consent of BioPro, Inc.



BioPro, Inc.
2929 Lapeer Road, Port Huron, MI 48060, USA
info@bioproimplants.com | 1-810-982-7777
www.bioproimplants.com

