


# Nitinol Staple Fixation

35 staple options for a variety of hand & wrist procedures.

The patented BioPro® Clover Staple™ uses strategically placed bends to compress the legs and top dimension towards the center of the fusion site.

The BioPro® Memory Staple features an S-bridge design to provide uniform compression and allow for fast insertion.



## Applications




Metacarpal/Capitate Fusion  
Memory Staple

Carpometacarpal Fusion  
Memory Staple

Radioscapholunate Fusion  
Clover Staple™

Four Corner Fusion  
Clover Staple™

Scaphotrapeziotrapezoid Fusion  
Clover Staple™



- 3 and 4-leg designs specially designed for intercarpal fusions
- Centralized compression to help speed the fusion process
- 1mm profile requires minimal bone removal and reduces risk of radial impingement



- 17 standard sizes with leg lengths as short as 5mm
- S-bridge design ensures even compression across the fusion site
- Easy to use instrumentation for fast and direct visual insertion

## Ordering Information

Item #	System	Top Width	Leg Length	Wire Size
18259	4-Leg	10x8mm	10mm	1.5
18260	4-Leg	10x8mm	12.5mm	1.5
18261	4-Leg	10x8mm	15mm	1.5
18062	4-Leg	14x12mm	10mm	1.5
18063	4-Leg	14x12mm	12.5mm	1.5
18064	4-Leg	14x12mm	15mm	1.5
18065	4-Leg	18x16mm	10mm	1.5
18066	4-Leg	18x16mm	12.5mm	1.5
18067	4-Leg	18x16mm	15mm	1.5
18262	3-Leg	10x8mm	10mm	1.5
18263	3-Leg	10x8mm	12.5mm	1.5
18264	3-Leg	10x8mm	15mm	1.5
18068	3-Leg	14x12mm	10mm	1.5
18069	3-Leg	14x12mm	12.5mm	1.5
18070	3-Leg	14x12mm	15mm	1.5
18071	3-Leg	18x16mm	10mm	1.5
18072	3-Leg	18x16mm	12.5mm	1.5
18073	3-Leg	18x16mm	15mm	1.5

Item #	Bridge Width	Leg Length	Wire Size
17637	7mm	5mm	1.2mm
17628	7mm	5mm	1.5mm
17638	9mm	7mm	1.2mm
17629	9mm	7mm	1.5mm
17630	11mm	8mm	1.5mm
17631	11mm	10mm	1.5mm
17634	11mm	15mm/13mm	1.5mm
17635	11mm	17mm/15mm	1.5mm
17636	11mm	19mm/17mm	1.5mm
17632	13mm	10mm	1.5mm
17633	15mm	12mm	1.5mm
17625	15mm	12mm	2.0mm
17626	18mm	18mm/15mm	2.0mm
17627	20mm	20mm	2.0mm
17622	20mm	20mm	2.0x3.0mm
17623	25mm	22mm	2.0x3.0mm
17624	30mm	30mm	2.0x3.0mm

This content is provided as an educational tool only and is not meant as medical advice in the usage of specific BioPro products. A healthcare professional must use their professional judgment in making any final determinations in product usage and technique. The product's Instructions for Use, should always be reviewed prior to surgery. Postoperative management is patient-specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level or outcomes. This information does not constitute medical, legal, or any other type of professional advice and should not be relied upon as such. It is not to be redistributed, duplicated, or disclosed without the express written consent of BioPro, Inc.



BioPro, Inc.  
2929 Lapeer Road, Port Huron, MI 48060, USA  
info@bioproimplants.com | 1-810-982-7777  
[www.bioproimplants.com](http://www.bioproimplants.com)



intertek

MKT131 01