



First MPJ Hemi



BRM Silktoe

# Hallux Rigidus Solutions

Our hallux rigidus portfolio supports every stage of arthritis, giving surgeons proven solutions from a single company. From joint-preserving procedures to definitive reconstruction, each option is designed to improve patient outcomes.



Accu-Cut® Guides



MTP Fusion Plates

## Stage 2/3



The BioPro® MPJ Hemi Implant is the only Hemi backed by 70+ years of clinical data. Our low-profile, press-fit implants are designed to replace the articular surface of the proximal phalanx.

- ✓ 20+ year survivorship data<sup>1</sup>
- ✓ +95% implant survivorship on average\*
- ✓ 97% of patients would recommend the procedure<sup>2</sup>
- ✓ Immediate weight-bearing
- ✓ ≤6 weeks return to activities<sup>2</sup>



The Accu-Cut® Youngswick and Chevron Osteotomy Guide system that allows for precise decompression osteotomies. The all-inclusive sterile packaged system includes everything needed to perform the osteotomy.

- ✓ 1-3mm decompression slots
- ✓ All-inclusive sterile packaged system
- ✓ Accurate decompression, every time.

## Stage 3/4



The BRM SilkToe® Implant is a new generation flexible silicone spacer for total first MTP joint arthroplasty.

- ✓ Maintains motion at later stage arthritis
- ✓ Left+Right orientation for anatomic reconstruction
- ✓ Latest generation silicone for durability.
- ✓ Immediate weight-bearing



The BioPro First MTP Fusion Plate is an anatomic plate with built in valgus offset and 5 or 10 degree dorsiflexion.

- ✓ Available in both standard and revision MTP plates
- ✓ Left+Right orientation for anatomic reconstruction
- ✓ Low profile plate with curvature to reduce irritation
- ✓ 2.8mm distal screws and 3.0/3.5 proximal screws



BioPro, Inc.  
2929 Lapeer Road, Port Huron, MI 48060, USA  
info@bioproimplants.com | 1-810-982-7777  
[www.bioproimplants.com](http://www.bioproimplants.com)

© 2026 BioPro, Inc. All Rights Reserved. MKT166 00

